

Connettori retti e a "Y"


Scopo

Per raccordare tubi (aspirazione, ecc.)

Specifiche

- Prodotto in polipropilene trasparente.
- Non invasivo.
- Non contiene lattice.
- Non sterile.

CE
CLASSE I
REGOLA 1



REF / ID	d1 mm	D2 mm	d3 mm	D4 mm	L1 mm	L2 mm	Confezione (pz.)
20075101 alias 8M130005K	4	6	4	6	47	-	Conf. 200 Box 2400
20075103 alias 8M130006K	7	10,5	7	10,5	55	-	Conf. 200 Box 2400
20075118 alias 8M100007K	7	14	7	14	90	-	Conf. 50 Box 750
20075105 alias 8M130007K	10	14	10	14	57	-	Conf. 200 Box 2000
20075107 alias 8M130002K	4	7	7	10	55	-	Conf. 200 Box 3000
20075121 alias 8M130012K	6	8	6	8	32	32	Conf. 100 Box 1500
20075123 alias 8M130013K	9	10	9	10	30	30	Conf. 100 Box 1500

Distribuito da:

Sim italia s.r.l.
SURGICAL INTERNATIONAL MARKETING

Società certificata
UNI CEI EN ISO 13485:2021
Via Gramsci, 9 - 40018 San Pietro in Casale (BO)
Tel. +39 051.817902 +39 051.817850
www.simitalia.it E-Mail: info@simitalia.it

95

Revisione del 22.01.2026